

Curadorm Max

The stable nursing bed for plus-size patients

Max.
patient
weight
350 kg

Safe
working load
up to 390 kg

Available in
4 widths and
3 lengths

Height
adjustment
range
40 - 80 cm



◀ The system cut-out at the head and foot end is ideal for storing bedding.

▶ The lying surface with metal straps is not only extremely robust, but can also be disinfected very quickly and easily.




The **Curadorm Max** offers maximum stability for daily care use in the bariatric/plus-size sector. The highly efficient lifting lever system and the lying surface with metal straps provide patients weighing up to 350 kg a safe bed with all the comforts, from the individually adjustable mattress platform elements to the hand control with individually lockable functions.

Product details

- | 4 widths: 90, 100, 120 or 140 cm
- | 3 lengths: 200, 210, or 220 cm
- | Under-bed clearance: 12 cm (height increase possible, please inquire)
- | 4-section lying surface with sturdy metal straps
- | Patient lifting pole with grab handle (optional)
- | Hand control with 4 functions and single function lock
- | Powerful Linak motor system
- | Individually braked heavy duty swivel castors

Product images



The product images on this page show the design of the standard nursing bed in the 120 cm variant. The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as side rail coverings made of imitation leather or fabric. In addition, the 90/100 cm, 120 cm and 140 cm versions of the bed differ visually in terms of their head/footboards, and in some cases in terms of the type of lying surface. If you have any questions regarding an order, please do not hesitate to [contact us](#) .

Technical data

- | Mattress width: 90 – 140 cm
- | Mattress length: 200 – 220 cm
- | Bed height: 40 – 80 cm
- | Max. patient weight: up to 350 kg
- | Safe working load: up to 390 kg

Extra equipment / optional

- | Width: 90 cm
- | Width: 100 cm
- | Width: 120 cm
- | Width: 140 cm
- | Length: 200 cm
- | Length: 210 cm
- | Length: 220 cm
- | Patient weight: 220 kg / Safe working load: 260 kg
- | Patient weight: 260 kg / Safe working load: 300 kg
- | Patient weight: 350 kg / Safe working load: 390 kg
- | **Removable headboard (NS-048-5)**
- | **Removable footboard (NS-048-5)**
- | **Side rail cover imitation leather (NS-049-1)**
- | **Side rail cover fabric (NS-049-0)**
- | **Side rail height extension**
- | **Infrared control**
- | **Patient lifting pole 3-legged with grab handle (NS-261-0)**
- | **Reinforced patient lifting pole with grab handle for higher patient weight (VS-520-0)**
- | **Patient lifting pole with grab handle**
provides secure and stable support to make it easier to sit up or change position in bed

ISKMed UK Limited · M +44 (0)7759 448 118
greg.whelan@iskomed.com · www.iskomed.co.uk



ISKO Koch GmbH · Egerländer Str. 28 · 95448 Bayreuth (GER)
Tel +49 921 150845-0 · info@iskomed.de · www.iskomed.de