

# Multidorm Flex

The standing bed for your home

Max.  
patient  
weight  
135 kg

Safe  
working load  
175 kg

Dimensions:  
90 x 200 cm

Height  
adjustment  
range  
40 - 80 cm



◀ The Multidorm Flex can be controlled via the handset or the control panel on the bed.

▶ If the patient is in the standing position in the event of a power failure, the bed can be lowered using the manual emergency release.





The **Multidorm Flex** is primarily used for passive standing training in home care and is often used for the early mobilisation of patients. The standing position serves to strengthen the circulation, improve breathing, prevent osteoporosis and drain the upper urinary tract, among other things.

## Product details

- | Single bed with standing function
- | Under-bed clearance: 11 cm
- | Standing function and TV position individually programmable
- | Patient lifting pole with grab handle incl. central receptacle (optional)
- | Integrated manual emergency lowering
- | Belt system with three high-quality support straps
- | Central braking system

## Product images



The product images on this page show the H-variant of the standing bed (with wooden side rail). The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as a wooden or plexiglass work and therapy table. The O-version (without side rails) of the bed is also visually very different due to modified wooden elements. If you have any questions regarding the equipment options, please do not hesitate to [contact us](#) .

## Technical data

- | Mattress width: 90 – 100 cm
- | Mattress length: 200 – 220 cm
- | Bed height: 40 – 80 cm
- | Max. patient weight: 135 kg
- | Safe working load: 175 kg

## Extra equipment / optional

- | Width: 90 cm
- | Width: 100 cm
- | Length: 200 cm
- | Length: 210 cm
- | Length: 220 cm
- | H-Variant (wooden side rails)
- | O-Variant (without side rails)
- | Metal side rail (for O-Variant)
- | Foot plate angle adjustable
- | Patient lifting pole with grab handle incl. central receptacle (SB-024-3)
- | Work and therapy table (wood)
- | Work and therapy table (plexiglass)
- | Armrests (SB-031-0)
- | Shock positioning (head down position) (SB-036-0)
- | Battery pack (SB-041-0)
- | Head fixation (SB-181-0)
- | Belt system made of synthetic leather (SB-032-5)
- | Removable headboard
- | Infrared control



**ISKomed UK Limited** · M +44 (0)7759 448 118  
greg.whelan@iskomed.com · [www.iskomed.co.uk](http://www.iskomed.co.uk)

**ISKO Koch GmbH** · Egerländer Str. 28 · 95448 Bayreuth (GER)  
Tel +49 921 150845-0 · info@iskomed.de · [www.iskomed.de](http://www.iskomed.de)