

Belladorm Comfort / Pro

The nursing bed for the cost-conscious

Max.
patient
weight
170 kg

Safe
working load
210 kg

Available in
2 widths
80 or 90 cm

Height
adjustment
range
20 - 60 cm or
40 - 80 cm



◀ The Belladorm Pro is supplied with wooden panelling for the head and foot area to cover the lifting motors.

▶ The functional stand-up aid was specially developed for the Belladorm series.



At ISKMed, the **Belladorm series** stands above all for high stability and easy access for patient lifting hoists and stand aids. With the lifting system integrated at the head and foot end, there are no obstacles under the bed to interfere when lifting patients in and out. The care bed is registered as an aid with the health and care insurance companies in Germany.

Product details


- | Single bed lying surface 80 - 90 x 200 cm
- | Available in 4 versions: Comfort without wood panelling, Pro with concealed lift motors and each Low version with lower entry
- | 4-section lying surface with metal grid support
- | Patient lifting pole with grab handle (optional)
- | Lower leg support with ratchet function - leg elevation possible via control by patient
- | Electric adjustment for height, head and leg section
- | Hand control with 5 functions and locking options
- | Trendelenburg and anti-Trendelenburg positioning possible
- | Powerful motor system (24V) with external power supply unit
- | Individually braked swivel castors



Product images



The product images on this page show different variants:

Image 1 shows the “comfort” variant with a height adjustment range of 40 - 80 cm. Image 2 shows the variant “pro” with a height adjustment range of 20 - 60 cm. On image 3 you can see the bed with the optional stand-up aid. The 4th image shows the deluxe wooden cases for both variants “comfort” and “pro”. The final product changes optically with optional equipment (“Extra equipment/optional” on next page), such as side rail coverings made of imitation leather or fabric. If you have any questions regarding an order, please do not hesitate to to [contact us](#) .

Technical Data

- | Mattress width: 80 – 90 cm
- | Mattress length: 200 – 220 cm
- | Bed height: 20 – 60 cm / 40 – 80 cm
- | Max. patient weight: 170 kg
- | Safe working load: 210 kg

Extra equipment / optional

- | Width: 80 cm
- | Width: 90 cm
- | Length: 200 cm
- | Length: 220 cm
- | Height adjustment range: 20 – 60 cm
- | Height adjustment range: 40 – 80 cm
- | Standard wood cladding
- | Deluxe wooden case
- | Easy EX light stand-up aid (NS-035-LC)
- | Side rail cover imitation leather (NS-049-1)
- | Side rail cover fabric (NS-049-0)
- | Side rail height extension
- | Infrared control
- | Patient lifting pole 3-legged with grab handle (NS-261-0)
- | Patient lifting pole with grab handle
provides secure and stable support to make it
easier to sit up or change position in bed



ISKMed UK Limited · M +44 (0)7759 448 118
greg.whelan@iskomed.com · www.iskomed.co.uk

ISKO Koch GmbH · Egerländer Str. 28 · 95448 Bayreuth (GER)
Tel +49 921 150845-0 · info@iskomed.de · www.iskomed.de