

Multidorm Mobil

The nursing bed with electric folding function

special beds



Max.
patient
weight
135 kg

Safe
working load
175 kg

Dimensions:
90 x 200 cm

Height
adjustment
range
40 – 80 cm



◀ After the automatic unfolding process, the Multidorm Mobil is a fully-fledged nursing bed that complies with the requirements of the ISO 60601-2-52 standard.



The **Multidorm Mobil** is perfect for quick deployment in the home care setting being folded and unfolded electrically at the touch of a button. When folded its compact size also offers ease of storage. Thanks to the patented side rail solution, no tools are required to set up and dismantle the bed.

Product details


- | Mobile single bed with 90 cm width
- | Under-bed clearance: 15 - 55 cm
- | Trendelenburg and anti-Trendelenburg positioning possible
- | Individually braked and very robust castor system
- | Particularly fast assembly and disassembly due to integrated folding motor
- | Patient lifting pole with grab handle (optional)
- | Integrated holding brackets for transporting the side rails
- | Patented side rail pick-up system
- | Easy to transport and store thanks to small dimensions
- | Single handed installation



Product images



special beds

The product images on this page show the standard version of the care bed. The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as a side rail cover made of imitation leather or fabric. If you have any questions regarding an order, please do not hesitate to [contact us](#) .

Technical data

- | Mattress width: 90 cm
- | Mattress length: 200 cm
- | Bed height: 40 – 80 cm
- | Max. patient weight: 135 kg
- | Safe working load: 175 kg

Extra equipment / optional

- | Side rail extension (NS-058-0)
- | Side rail cover imitation leather (NS-049-1)
- | Side rail cover fabric (NS-049-0)
- | Patient lifting pole with grab handle (NS-024-6)



ISKMed UK Limited · M +44 (0)7759 448 118
greg.whelan@iskomed.com · www.iskomed.co.uk

ISKO Koch GmbH · Egerländer Str. 28 · 95448 Bayreuth (GER)
Tel +49 921 150845-0 · info@iskomed.de · www.iskomed.de