

# Curadorm Pro

The nursing bed for the ward

Max.  
patient  
weight  
170 kg

Safe  
working load  
210 kg

Available in  
3 lengths:  
200, 210 and  
220 cm

Height  
adjustment  
range  
32 - 82 cm



nursing beds



◀ The central brake system locks all four casters simultaneously, so that the bed can be easily and securely positioned directly against the wall.

▶ Adjustable mattress holders ensure high flexibility and additional safety when securing different support systems.




The **Curadorm Pro** is a nursing bed with scissor lift system primarily for professional use in nursing homes and rehabilitation centres. The particularly large adjustment range and the high-quality castor system with central brake are especially appreciated by carers in everyday nursing care.

## Product details

- | Single nursing bed in 80 or 90 cm width
- | Under-bed clearance: 13 cm
- | 4-section, stable lying surface with hygienic metal grid
- | Patient lifting pole with grab handle (optional)
- | Hand control with 4 functions and single function lock
- | System cut-out on headboard and footboard for storage (e.g. duvets)
- | High-quality Linak motor system
- | Clinic double castor system with central brake

## Product images



The product images on this page show the standard version of the nursing bed. The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as plastic or leather side rail coverings. If you have any questions regarding an order, please do not hesitate to [contact us](#) .

## Technical data

- | Mattress width: 80 – 90 cm
- | Mattress length: 200 – 220 cm
- | Bed height: 32 – 82 cm
- | Max. patient weight: 170 kg
- | Safe working load: 210 kg

## Extra equipment / optional

- | Width: 80 cm
- | Width: 90 cm
- | Length: 200 cm
- | Length: 210 cm
- | Length: 220 cm
- | Easy EX stand-up aid (NS-035-0)
- | Easy EX light stand-up aid (NS-035-LC)
- | Removable headboard (NS-048-5)
- | Removable footboard (NS-048-5)
- | Side rail cover imitation leather (NS-049-1)
- | Side rail cover fabric (NS-049-0)
- | Side rail height extension
- | Infrared control
- | Patient lifting pole 3-legged with grab handle (NS-261-0)
- | Patient lifting pole with grab handle  
provides secure and stable support to make it  
easier to sit up or change position in bed



**ISKMed UK Limited** · M +44 (0)7759 448 118  
greg.whelan@iskomed.com · [www.iskomed.co.uk](http://www.iskomed.co.uk)

**ISKO Koch GmbH** · Egerländer Str. 28 · 95448 Bayreuth (GER)  
Tel +49 921 150845-0 · info@iskomed.de · [www.iskomed.de](http://www.iskomed.de)